KEEP IN TOUCH
TELEHEALTH SOLUTIONS

GuideMe

**Manufacturer:**

AIT Austrian Institute of Technology GmbH Giefinggasse 4
1210 Wien, Austria
Date of manufacture: 2025

KITGuideMe is a medical device (class IIa).

**Serial number:**

KITGuideMe 1.0.1

MDC MEDICAL DEVICE CERTIFICATION GMBH
Kriegerstrasse 6, 70191 Stuttgart, Germany

**Use-by date:**

2027-10

**UDI-DI**

(01)09120109830252

UDI-PI

Software Version (8012) 1.0.1

Caution: Please read paragraphs marked with this symbol carefully as they contain important safety information.

Document Date: 2025-04-28

Document Version: 1.1

en

INSTRUCTIONS FOR USE

KITGUIDEME VERSION 1.0.1

2

de

GEBRAUCHSINFORMATION

KITGUIDEME VERSION 1.0.1

13

Instructions for Use for Health Professionals

Table of contents

1 Principles of Operation.....	3
1.1 Available Features	3
2 Intended Use	3
3 Intended User Group.....	4
4 Use Environment	4
5 Lifespan.....	5
6 Indication	5
7 Contra Indication	5
8 Patient Population	5
9 Residual Risks and Side-Effects.....	5
10 Limitations and Exclusions	5
11 Calculation of the KITGuideMe Result.....	6
12 Structure of the KITGuideMe result	6
13 Interpretation of the KITGuideMe result	7
13.1 Visualisation of the medication status in relation to the predefined target dose	7
13.2 Detailed information about drug class	8
14 Drug classes and active components considered with KITGuideMe.....	8
15 Reasons to deviate from target value	10
16 Additional help	11
17 Commissioning, maintenance and repair measures	11
18 Application environments for KITGuideMe.....	11
19 Report feedback and problems	11
20 Printed version of this instruction for use.....	11
21 Explanation of the symbols in the instructions for use	12

1 Principles of Operation

The medical product KITGuideMe serves as an independent software service, which can be used by users as an add-on module for a data management system (DMS), in order to assist health professionals to evaluate the current medication status based on an underlying guideline. This can be achieved by comparing the prescribed daily dose of active components within a predefined drug class in comparison to guideline-recommended or custom daily target doses using predefined rules. The daily results of active components are summarized in their dedicated drug classes and their overall drug-class status is indicated in percent and with a status-dependent indicator. To enable the deviation of the prescription to the guideline-recommendations, health professionals can state predefined reasons (e.g. contraindications) derived from the guidelines. KITGuideMe receives the submitted data through an API from the DMS, analyses the data according to the predefined rules and returns the results back to the DMS. The results are provided as HTML or in a machine readable format to provide DMS specific visualisation (look and feel). KITGuideMe does not store any data.

1.1 Available Features

KITGuideMe is able to store, interpret and apply different guidelines. Each of these possible guidelines within the software product is called "feature". KITGuideMe has following features available in this version:

Feature	Description
"2021 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure" and its update "2023 Focused Update of the 2021 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure".	The current prescription status of the daily dose of a prescribed drug in relation to a predefined daily target dose is visualized for the professional user. Active components, drug classes and reasons to deviate from a guideline were derived from the "2021 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure" and its update "2023 Focused Update of the 2021 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure". For the calculation, only the total daily dose is considered, not the recommended intake schema by the guideline.

2 Intended Use

KITGuideMe is intended to assist health professionals to evaluate the current medication status of the prescribed daily dose of active components, in particular RAASi, Beta-Blocker, MRA and SGLT2-Inhibitors, within a predefined drug class in comparison to predefined daily target doses using predefined rules. KITGuideMe is a clinical decision support tool, that supports the health professional to optimize the prescribed dose by comparing the medication status of a patient against predefined target values. These target values can either be derived from an existing guideline (e.g. ESC Guidelines for Chronic and Acute Heart Failure) or optional customized target values. Those customized values have to be predefined by a qualified health professional (e.g. physician).

The health professional can deviate from the recommended doses. These deviations are exclusively explained by predefined reasons (e.g. contraindications) and are only applicable for predefined drug classes. These predefined reasons are evidence-based recommendations of a guideline.

KITGuideMe can match the prescribed medications and active components into the following drug

classes:

- Renin-angiotensin-aldosterone-system inhibitors (RAASi)
- Beta-Blocker,
- Mineralcorticoid receptor anatagonists (MRA) or
- Sodium-glucose cotransporter-2 (SGLT2)- inhibitors

The medical device does not give any diagnoses, diagnostic recommendations or prognoses about the condition of the patient. It remains the health professional's responsibility to decide on diagnoses and treatment. The recommended target doses do not override the individual responsibility of health professionals to make appropriate and accurate decisions. These decisions shall be made in consideration of each patient's health condition and in consultation with that patient or the patient's caregiver where appropriate and/or necessary. It is also the health professional's responsibility to verify the rules and regulations applicable in each country to drugs and devices at the time of prescription.

The software module does not give any diagnoses, diagnostic recommendations or prognoses about the condition of the patient. It remains the physician's responsibility to decide on diagnoses and treatment.

The software is not intended to control the interactions of the prescribed drugs.

3 Intended User Group

The intended users are health professionals like physicians or nursing staff with a relevant medical background. The intended user must have the physical and mental ability to handle a PC and its standard functions. No special user training is required for the safe use of KITGuideMe. No gender is excluded.

The interpretation of the results and derivation of therapeutic measures is intended for qualified health professionals only.

4 Use Environment

KITGuideMe is operated as an independent software service, which must be accessible over a network for the DMS. KITGuideMe is operated on a server (physical or virtual).

The interaction with KITGuideMe takes place in a broader context, more precisely the sighting of the results, both in medically used areas, e.g. hospitals, inpatient or outpatient therapeutic facilities, health resorts or rehabilitation institutions, as well as rooms used for non-medical purposes, e.g. by mobile care staff or by practising physicians during house calls.

Please make sure to use the software only in a trustworthy environment. The media owner is responsible for the compliance with requirements of the general data protection regulation (EU).

The software module is not intended for emergency medical care or intensive treatment.

5 Lifespan

In principle software has no limitation on the lifespan. However, the software must not be used with operating system and the supported Python Version that is no longer maintained. Therefore, the lifespan of the software version ends on support time of the operating system and the supported Python Version. The lifespan for this version of the product ends on 2027-10-31 (October 31st 2027).

As of this date, no more error corrections will be performed for the software and its version may no longer be used. The service life of subsequent versions may be defined differently.

6 Indication

KITGuideMe can be used wherever a dose of active components, within specific drug classes, in particular RAASi, Beta-Blocker, MRA and SGLT2-Inhibitors, are prescribed to a patient (e.g. HFrEF- Heart Failure with reduced ejection fraction).

The software module does not alert the user if the predefined daily dose is not reached.

7 Contra Indication

The software module is not intended for emergency medical care or intensive treatment.

8 Patient Population

The predefined daily doses are defined for the entire patient cohort of a particular disease management program or a subset of patients with a distinct medical condition. There are no restrictions regarding gender, racial and ethnic origin.

9 Residual Risks and Side-Effects

No known residual risks or undesirable side-effects.

10 Limitations and Exclusions

The following limitations and exclusions apply for KITGuideMe:

- The software module does not give any diagnoses, diagnostic recommendations or prognoses about the condition of the patient. It remains the physician's responsibility to decide on diagnoses and treatment.
- The software module is not intended for emergency medical care or intensive treatment.
- The software module does not alert the user if the predefined daily dose is not reached.
- The software is not intended to control the interactions of the prescribed drugs.

- Users that are physically or mentally not fit to use the software module.
- The interpretation of the results and derivation of therapeutic measures is intended for qualified health professionals only.

11 Calculation of the KITGuideMe Result

KITGuideMe receives medication data in a structured format from a DMS. It validates the daily intake or intake schema, active components, dose of active components, unit of active components, reasons for deviating from a guideline and custom target values if they are set. These validations may lead to errors which are described in the result.

If custom target values are not plausible, the software will use the guideline-recommended value as a fallback for calculation. The phrase “fallback” is then also shown in the HTML report.

The {daily relative dose of active components} is calculated as the {daily prescribed dose of an active component} divided by the {daily target dose} and multiplied by 100 as the result is displayed in percent.

The {daily prescribed dose of an active component} is calculated as the {dose of an active component} multiplied by the {sum of the intake schema of prescribed doses}.

Dosages can be provided in the units of [µg], [mg] and [g], however, they are transformed to [mg] for calculation and the software does only return values in [mg].

If more than one drug or active component per drug class is prescribed, the individual daily relative doses are summed up. However, if more than one drug or active component per drug class is prescribed, a hint will be provided to draw attention to this fact.

Relative dose values are rounded to integer within the human readable HTML report and to two decimal places in the machine-readable JSON format.

12 Structure of the KITGuideMe result

KITGuideMe presents the analysis results in a report that can be opened and viewed on the PC with a Web Browser. In addition, KITGuideMe offers an interface for the corresponding data management system, which makes the analysis results available electronically in structured form for visualization in the DMS. This function is not part of this instructions for use and is described in the integration manual of KITGuideMe for commissioning by an IT technician.

Please note that the visualisation may slightly deviate depending on the data management system in which KITGuideMe is integrated.

The interpretation of the results and derivation of therapeutic measures is intended for qualified health professionals only.

Following information is provided in the result report:

- Medication status in relation to the predefined target dose
 - Visualisation of a medication status for different drug classes in form of a bar chart.
- Detailed information about drug class
 - Further information about the drugs and active components within a drug class compared to the target values.
- General information

- This section contains a short description, the indicators, colors. Furthermore, the active components relevant for each drug class, as well as possible reasons for each drug class - based on the guideline - are provided.

13 Interpretation of the KITGuideMe result

13.1 Visualisation of the medication status in relation to the predefined target dose

KITGuideMe visualizes a medication status for different drug classes in form of a bar chart.

Each bar shows the prescribed medication in relation to the recommended target dose per patient per day. A relative dose of 100% would mean that the prescribed dose is equal to the guideline-recommended target dose for a drug class per day.

For each drug class, the bar is coded with one of three colors:

Green: Relative dose of the drug class at 100% of the target dose or a reason for deviating from the target dose has been given.

Yellow: Up-titrate of a patient's drug class (if "Up-titrate patient's drug class" was set as the only reason).

Red: Relative dose less than or greater than 100% and no reason has been given.



The displayed colors do not depict a recommended action but should be interpreted as a form of visualization of different status.

Furthermore, each bar has an indicator next to the percentage, they mean the following:

- OK Ok: Guideline-recommended or custom target dose value(s).
- NOK Not Ok: Dose different to guideline-recommended or custom target dose value(s).
- R Reason: reasons given for deviation from the guideline.
The reason will be displayed in the detail section of the result.
- T Up-titration: Indicates the imminent up-titration of an active component.
- H Hint: Problem caused by the entered input data by a professional user
The hint will be displayed in the detail section of the result
- E Error: Problem caused by the software
The error will be displayed in the detail section of the result

In the case of combined active components (e.g. Sacubitril/Valsartan), the relative dose of drugs with combined active components is calculated by the mean value of the individual relative dose values of each of the contained active components (within the drug).

Please note that the visualisation may slightly deviate depending on the data management system in which it is integrated.

Example visualisation HTML:**Example visualisation within a DMS:****13.2 Detailed information about drug class**

For each drug class a detailed Information is available. This contains:

- Daily relative Dose
- Optional indicators with details like: provided reasons, hints or errors
- A table with:
 - The relevant drug for the drug class
 - The active component relevant for the drug class
 - The absolute dose of the active component in milligram [mg]
 - The source of the target dose (guideline or custom)
 - The absolute target dose based on the given source in [mg]
 - The calculated relative dose

In case of multiple drugs per drug class, the the relative dose of each active component is summed up to a total relative dose within the drug class. In this special case the bar is displayed in red and a hint is displayed.

14 Drug classes and active components considered with KITGuideMe

This version of KITGuideMe considers drugs and the relevant active components due to the "2021 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure" and its update "2023 Focused Update of the 2021 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure" from following drug classes:

- Renin-angiotensin-aldosterone-system inhibitors (RAASi)

- Beta-Blocker
- Mineralcorticoid receptor antagonist (MRA)
- Sodium-glucose cotransporter-2 (SGLT2)-inhibitors

The active components are:

Active Component	Drug class	Drug sub class	Daily Target [mg]	Custom Min [mg]	Custom Max [mg]
Captopril	RAASi	ACE-I	150	18.75	150
Enalapril	RAASi	ACE-I	40	5	40
Lisinopril	RAASi	ACE-I	35	2.5	35
Ramipril	RAASi	ACE-I	10	5	10
Trandolapril	RAASi	ACE-I	4	0.5	4
Candesartan	RAASi	ARB	32	4	32
Losartan	RAASi	ARB	150	50	150
Valsartan	RAASi	ARB	320	80	320
Sacubitril/valsartan	RAASi	ARNI	194/206	48/52	195/206
Bisoprolol	Beta-blockers	Beta-blockers	10	1.25	10
Carvedilol	Beta-blockers	Beta-blockers	50	6.25	100
Metoprolol	Beta-blockers	Beta-blockers	200	12.5	200
Nebivolol	Beta-blockers	Beta-blockers	10	1.25	10
Eplerenone	MRA	MRA	50	25	50
Spironolactone	MRA	MRA	50	12.5	50
Dapagliflozin	SGLT2i	SGLT2i	10	10	10
Empagliflozin	SGLT2i	SGLT2i	10	10	10

Additional information on the table of active components above:

The Daily Target, Custom Min and Custom Max values are derived from the "2021 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure" and its update "2023 Focused Update of the 2021 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure".

Custom Min defines the minimal value allowed for custom target dose values and is derived from the starting dose of the respective active component (of the guidelines). Custom Max defines the maximum value allowed for custom target dose values and is derived from the target dose of the respective active component (of the guidelines).

In case of target dose value-ranges in the guidelines, the higher value was used for Daily Target and Custom Max values.

In case of starting dose value-ranges in the guidelines, the lower value was used for the Custom Min value.

Sacubitril/valsartan: Custom Min: may have an optional lower starting dose of 24/26 mg b.i.d. for those with a history of symptomatic hypotension.

Carvedilol: Custom Max: A maximum dose of 50 mg twice daily can be administered to patients weighing over 85 kg.

Spironolactone: Custom Min: optional starting dose of 12.5 mg in patients where renal status or hyperkalaemia warrant caution.

Custom target dose values and ranges only support the unit mg.

The custom_max [mg] of Sacubitril/valsartan was altered, so that it is higher than the recommended target dose (195/206 vs. 194/206) due to unreachable target dose values of 194.4 / 205.6 (mg) when using the medication

database of the Austrian pharmacological society.

15 Reasons to deviate from target value

Following reasons to deviate from a target value can be given at KITGuideMe. They are taken from the "2021 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure" and its update "2023 Focused Update of the 2021 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure".

Drug class	Drug sub class	Reason
RAASi	ACE-I, ARB, ARNI	History of angioedema
RAASi	ACE-I, ARB, ARNI	Known bilateral renal artery stenosis
RAASi	ACE-I, ARB, ARNI	Pregnancy/risk of pregnancy
RAASi	ARNI	Breastfeeding period
RAASi	ACE-I, ARB, ARNI	Known allergic reactions/other adverse reactions (drug-specific)
RAASi	ACE-I, ARB, ARNI	eGFR 30 mL/min/1.73 m ²
RAASi	ACE-I, ARB	Significant renal dysfunction (creatinine > 2.5 mg/dL)
RAASi	ACE-I, ARB, ARNI	Symptoms of hypotension or a SBP <90 mmHg
RAASi	ACE-I, ARB, ARNI	Significant hyperkalaemia (> 5.0 mmol/L)
RAASi	ACE-I, ARB, ARNI	Up-titration
Beta-blockers	Beta-blockers	Known allergic reactions/other adverse reactions (drug-specific)
Beta-blockers	Beta-blockers	Second- or third-degree AV block (in the absence of a permanent pacemaker)
Beta-blockers	Beta-blockers	Critical limb ischaemia
Beta-blockers	Beta-blockers	Asthma (not COPD)
Beta-blockers	Beta-blockers	Severe exacerbation of heart failure
Beta-blockers	Beta-blockers	Bradycardia (< 50 b.p.m.)
Beta-blockers	Beta-blockers	Hypotension (SBP < 90 mmHg)
Beta-blockers	Beta-blockers	Up-titration
MRA	MRA	Known allergic reactions/other adverse reactions (drug-specific)
MRA	MRA	Significant hyperkalaemia (> 5.0 mmol/L)
MRA	MRA	Significant renal dysfunction (creatinine > 2.5 mg/dL)
MRA	MRA	eGFR <30 mL/min/1.73 m ²
MRA	MRA	Up-titration
SGLT2i	SGLT2i	Known allergic reactions/other adverse reactions (drug-specific)
SGLT2i	SGLT2i	Pregnancy/risk of pregnancy and breastfeeding period
SGLT2i	SGLT2i	eGFR 20 mL/min/1.73 m ²
SGLT2i	SGLT2i	Symptoms of hypotension or SBP < 95 mmHg
SGLT2i	SGLT2i	Acute (genito-urinary) infection

16 Additional help

If you need additional help, please contact the manufacturer via the e-mail address kitguideme@ait.ac.at.

17 Commissioning, maintenance and repair measures

KITGuideMe requires no special maintenance. However, the manufacturer may issue updates at irregular intervals. The manufacturer or distributor will inform the user about these updates. It is recommended that updates from the web browser manufacturer are installed promptly. The user must ensure that the web browsers used meet the minimum requirements defined below: Microsoft Edge (at least version 122), Firefox (at least version 123), Chrome (at least version 124).

The requirements for KITGuideMe concerning the installation environment of the system operator are listed in a dedicated document for the operator (KITGuideMe – Installation Guide and operation Manual).

18 Application environments for KITGuideMe

KITGuideMe can only be used in combination with a DMS (Data Management System) that provides KITGuideMe with medication data of the patient and possible reasons. The results can be displayed directly in the DMS via a technical interface or in form of a HTML report. For further information please contact kitguideme@ait.ac.at.

19 Report feedback and problems

Please report feedback and problems with KITGuideMe to the distributor where the software was purchased or directly to the manufacturer using the email address kitguideme@ait.ac.at.

Any serious incidents that occur in relation to the device should be reported directly to the manufacturer using the email address kitguideme@ait.ac.at and the competent authority of the member state in which the user and/or patient is established.

20 Printed version of this instruction for use

A printed version of the instructions for use can be requested directly from the manufacturer via the e-mail address kitguideme@ait.ac.at.

21 Explanation of the symbols in the instructions for use

Symbol	Meaning	Description
	Caution	Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.
	Manufacturer	This symbol refers to the manufacturer's address.
	Medical device	This symbol indicates that this product is a medical device.
	Serial number	This symbol indicates the serial number of the software to which this instruction manual applies.
	CE mark	This symbol indicates conformity with Regulation (EU) 2017/745 with the number of the Notified Body (0483- mdc medical device certification GmbH, Kriegerstrasse 6, 70191 Stuttgart, Germany).
	Use-by date	This symbol indicates the date after which the medical device may no longer be used.
	Unique Device Identifier	Indicates a carrier that contains unique device identifier information.

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